Description

Tetric EvoCeram® Bulk Fill is a light-curing, radiopaque nano-hybrid composite for direct restorations in posterior teeth.

Tetric EvoCeram Bulk Fill cures with light in the wavelength range of 400–500 nm and can be applied in layers of up to 4 mm.

Shades

Tetric EvoCeram Bulk Fill is available in the following 3 shades:

- Universal shades: ^{IV}A, ^{IV}B, ^{IV}W

Composition

The monomer matrix is composed of dimethacrylates (20-21 wt.%).

The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (79-81 wt.%). Additives, initiators, stabilizers and pigments are additional ingredients (<1.0 wt.%). The total content of inorganic fillers is 76–77% weight or 53–54% volume. The particle sizes of the inorganic fillers range between 40 nm and 3 µm.

Indication

- Restoration of deciduous teeth
- Restorations in the posterior region (Classes I and II, including the replacement of individual cusps)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Reconstructive build-up
- Extended fissure sealing in molars and premolars

Contraindication

The placement of Tetric EvoCeram Bulk Fill restorations is contraindicated:

- If a dry working field cannot be established, or the stipulated working technique cannot be applied.
- If a patient is known to be allergic to any of the ingredients of Tetric EvoCeram Bulk Fill.

Side effects

In individual cases, components of Tetric EvoCeram Bulk Fill may lead to sensitization. Tetric EvoCeram Bulk Fill should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp and cover it with a suitable cavity liner).

Interactions

Phenolic substances such as materials containing eugenol/clove oil inhibit the polymerization of materials based on methacrylate. Consequently, the application of such materials together with Tetric EvoCeram Bulk Fill must be avoided. Cationic mouth rinses and plaque disclosing agents as well as chlorhexidine may cause discolouration.

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Application

1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel enamel margins of anterior teeth. In the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes using rubber cups or rotary brushes. Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Selectively coat only very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently cover with a pressure-resistant cement (e.g. a glass-ionomer cement, such as Vivaglass® Liner). Do not cover the other cavity walls so that they can be used to support the bond with an enamel/dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® Universal, a single-component, light-cured adhesive for direct and indirect bonding procedures that features compatibility with all etching techniques.

7. Application of Tetric EvoCeram Bulk Fill

- For an optimum result, apply Tetric Evo Ceram Bulk Fill in layers of max. 4 mm and adapt with a suitable instrument (e.g. OptraSculpt).
- Sufficient exposure to the curing light prevents incomplete polymerization. For the recommendations regarding exposure time and light intensity see table 1.
- When using a metal matrix, additionally polymerize the composite material from the buccal or lingual/palatal aspect after removing the matrix, if no Bluephase polymerization light is used or the light probe cannot be ideally positioned, e.g. distance to the composite or diverging scattering angle.
- A flowable composite (e.g. Tetric EvoFlow® or Tetric EvoFlow® Bulk Fill) can optionally be used. This optional layer must be cured separately (please refer to the respective Instructions for Use).

8. Finishing / Checking the occlusion / Polishing

Remove excess material with diamond or tungsten carbide finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric EvoCeram Bulk Fill can be used in combination with Tetric EvoFlow Bulk Fill, Tetric EvoCeram and Tetric EvoFlow.
2. In the case of repairs, additional Tetric EvoCeram Bulk Fill can be directly applied to polymerized material. If the Tetric EvoCeram Bulk Fill restoration has already been polished, it must first be roughened and wetted with Adhese Universal before fresh Tetric EvoCeram Bulk Fill is applied.
3. Tetric EvoCeram Bulk Fill should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If Tetric EvoCeram Bulk Fill is applied from the Cavifil directly into the oral cavity of the patient, the Cavifil or the application cannula of the syringe should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
5. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.
6. The recommended increment thickness is based on hardness profile measurements.

Warning

Avoid contact of unpolymerized Tetric EvoCeram Bulk Fill with skin, mucous membrane and eyes. Unpolymerized Tetric EvoCeram Bulk Fill may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Do not use Tetric EvoCeram Bulk Fill after the indicated date of expiration.
- Date of expiration: see information on Cavifils, syringes and packages.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry and should be used strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.