Virtual[®]

EN Instructions for Use

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English

Description

Virtual impression materials are addition-reaction silicones (vinylpolysiloxanes) used to create fine detail impressions of dentition. Virtual impression materials are available in a variety of viscosities allowing dental professionals to choose the material and technique best suited for each individual case.

Colors

See table "Technical Data"

Material	Speed	Total Work Time [min:s]	Mouth Set Time* [min:s]
Extra Light Body [wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Light Body [wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Medium Body [wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Monophase [tray/wash]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Heavy Body [tray]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30
Putty [tray]	Regular Set	1:30	4:30
	Fast Set	1:00	2:30

* Minimum time the impression material should remain in the mouth before removal.

The processing times are reduced or prolonged depending on the prevalent room temperature.

Technical Data

	Extra Light Body	Light Body	Medium Body	Monophase	Heavy Body	Putty
Color	Beige	Beige	Beige	Blue	Blue	Blue
Classification ISO 4823 / ADA Spec. No. 19	Type 3 light-bodied consistency	Type 3 light-bodied consistency	Type 2 medium-bodied consistency	Type 2 medium-bodied consistency	Type 1 heavy-bodied consistency	Type 0 putty consistency
Mixing Ratio [Base:Catalyst]	1:1	1:1	1:1	1:1	1:1	1:1

Composition

Virtual impression materials are addition-reaction silicones containing vinylpolysiloxane, methylhydrogensiloxane, an organoplatinic complex, silica and food dyes.

Indication

The Virtual line of addition-silicone (vinylpolysiloxane) impression materials is recommended for use to create highly detailed impressions of the hard and soft tissues of the oral cavity.

- Final impression used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix from "wax-ups" or for treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Each viscosity provides specific attributes to meet the needs of the dental professional.

Contraindication

The use of Virtual impression materials is contraindicated if the patient is known to be allergic to any of the ingredients in Virtual impression materials.

Side Effects

None known to date

Interactions

The setting of vinylpolysiloxanes is inhibited by latex gloves.

Do not touch preparations or retraction cords with latex gloves. It is recommended that operators wash their hands thoroughly or use vinyl gloves in order to eliminate all traces of impurities, specifically when hand mixing putty. Other materials which may inhibit the set of impression materials include rubber dam, retraction cords and retraction cord fluids. If the operator suspects that the preparation has been contaminated, it is recommended that the preparation be rinsed and dried to eliminate all traces of impurities.

Application

– Putty (Regular and Fast Set)

Very high viscosity vinylpolysiloxane that provides easy mixing and the hydraulic forces necessary to extend wash material into the sulcus providing a detailed impression of prepared margins. The putty viscosity maintains outstanding detail reproduction when used alone. Recommended impression techniques include Putty/Wash One-Step and Putty/Wash Two-Step.

Important: Avoid wearing latex gloves when handling putty, putty jars, and dosing spoons (see note under interactions).

1. Take equal amounts of base (blue) and catalyst (white) using the color coded dosing spoons.
2. Hand mix the identical portions of Virtual Putty base and catalyst until you obtain an evenly colored mixture (approximately 30 seconds of mixing time). A proper mixture should have no streaks. Note that equal proportions of base and catalyst must be used in order to obtain proper setting times. A larger quantity of catalyst will not accelerate the setting

time. However, the setting time and curing depend on the temperature of the fingers/hand and the mixing intensity.

3. Place mixture in impression tray. It is recommended that the impression trays are pre-coated with Virtual Tray Adhesive.

Make sure to re-seal the Virtual Putty jars immediately after use using the lid with the correct color. Moreover, use the spoon with the correct color to measure out the material. Failure to do so can cause discoloration on the surface of the white Putty Catalyst paste. However, this discoloration does not impair the functions of the product.

- **Extra-Light Body, Light Body, Medium Body, Heavy Body, Monophase (Regular and Fast) Extra Light-Body (Beige):** Extra low viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Light-Body (Beige): Low viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Medium Body (Beige): Medium viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Heavy Body (Blue): High viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray material when using the Tray/Wash impression technique.

Monophase (Blue): Medium viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray and wash material when using the double-arch, single-phase impression technique.

- **Instructions for Cartridge Materials**

Loading the Cartridge Dispenser

(These instructions are written for use with the Virtual Manual Dispenser.)

See separate instructions for use if using the Virtual Automatic Dispenser.)

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

1. Press the black release lever located below the plunger on the back of the dispenser and pull plunger as far back as possible. (Fig. 1)
2. Lift the cartridge lock and insert the cartridge with the "V" shape on the cartridge base turned downwards. Lower the cartridge lock. (Fig. 2)
3. Once the cartridge is secure in the dispenser, remove the cartridge cap by turning it a 1/4 turn counter clockwise. Discard the cap. (Fig. 3)

4. It is VERY IMPORTANT to clear or bleed the cartridge PRIOR to applying the mixing tip. This will assure a proper mix and set time. Gently press on the dispenser lever until both components (base and catalyst) begin to flow out of the cartridge evenly onto a mixing pad.
5. Select the appropriate mixing tip:
 - a. The large mixing tips (BLUE colored base) are recommended for Heavy Body and Monophase impression materials
 - b. The small mixing tips (YELLOW colored base) are recommended for Extra-Light Body, Light Body and Medium Body impression materials
(The color of the mixing tip base and the color of the cartridge cap are identical)Insert the mixing tips in a clean motion so that the "V" shaped notch on the colored base is aligned with the "V" shaped notch on the housing of the cartridge. While gripping the colored base, and not the mixing tip, turn the mixing tip base 1/4 turn clockwise. (Fig. 4)
6. If necessary apply an intra-oral tip to the mixing tip (Fig. 5):
 - a. The large intra-oral tips (CLEAR) are to be used with the large mixing tips (BLUE base)
 - b. The small intra-oral tips (YELLOW) are to be used with the small mixing tips (YELLOW base)
7. Begin mixing by pressing gently on the lever of the dispenser.

Removing the cartridge

1. Lift the release lever to pull slide back. Lift the cartridge lock and remove cartridge.
2. Leave the mixing tip on the cartridge. The material will polymerize within the mixing tip, functioning as a natural cap.

– Pre-Treatment of Impression Trays (Tray Adhesives)

It is highly recommended that a tray adhesive (eg: Virtual Tray Adhesive) be used to reduce the chance of distortion when removing impressions from the mouth.

– Instructions for Virtual Tray Adhesive:

1. Be sure that all surfaces of the impression tray are oil-free, clean and dry.
Note: If the disinfection of the impression tray cannot be ensured, the tray adhesive must be filled into a dappen glass, for example, and applied using a disposable brush.
2. Apply a thin layer of Virtual Tray Adhesive using the brush provided to all surfaces of the impression tray (metal or plastic) that will come into contact with the impression material.
3. Allow the tray adhesive to dry for approximately 3 minutes (or blow dry with compressed oil-free air)
4. Replace bottle cap immediately.

– Cleaning of tray

Residue of adhesive material can be removed with alcohol.

– Note the warnings in the Virtual Tray Adhesive Instructions for Use!

– **Disinfection of impression**

Impressions made with Virtual impression materials can be immersed in a disinfection solution (glutaraldehyde 0.5%, benzalkonium chloride 0.5%). (Please observe the instructions of the manufacturer.) The disinfection does not affect the surface or the dimension.

– **Pouring Up Models**

The impression may be poured immediately after disinfection, or up to two weeks later, provided that the impression is stored at room temperature. Dimensional stability is guaranteed for 14 days, if stored appropriately. Virtual impression materials are compatible with all popular dental plasters on the market, e.g. Type 3: Elite® Model (Zhermack), Type 4: Fujirock® (G.C. International)

– **Galvanization**

Virtual impression can be silver or copper plated in a galvanic bath.

Special Notes

Virtual impression materials should be at room temperature (23 °C / 73.4 °F) when processed. Lower temperatures, e.g. when stored in the refrigerator, will prolong the time the impression material should remain in the mouth, while higher temperatures will reduce the processing time and the time the impression should remain in the mouth.

Vinylpolysiloxanes are chemically resistant. Unpolymerized materials may stain clothing.

Warnings

If uncured materials come into contact with the eyes, rinse the eye with copious amounts of water. If irritation persists, seek medical attention. In case of contact with the skin, wash affected areas with plenty of water.

Storage

- Storage temperatures: 2 – 28 °C / 36 – 82 °F
- Shelf life: See expiration date on label and packaging
- Keep away from direct heat sources!

Keep out of the reach of children!

For use in dentistry only!

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.