Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

Table with 10 columns: Au, Pt, Pd, Ag, Cu, Zn, Ni, In, Fe, Co. Values: <1.0, 53.8, 34.9, 0.2, 0.1, 0.1, 0.1, 0.1, 0.1, 0.1

Composition

Instructions for Use, Gebrauchsinformation, Instruções de Uso, Instruções de Uso, Brüksanvisning, Brüksanvisning, Käyttöohjeet, Käyttöohjeet, Brüksanvisning, Brüksanvisning, Produktionsformate, Produktionsformate

White, palladium-based dental ceramic casting alloy, Type 4

Ivoclar Vivadent Worldwide

- Australia: Ivoclar Vivadent Pty. Ltd., 1-5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174...
France: Ivoclar Vivadent SAS, Rue de la République, 723 Centre Empresarial Tambore CEP 06460-110 Barueri - SP...
Germany: Ivoclar Vivadent GmbH, Dr. Adolf-Schneider-Str. 2, 73479 Ellwangen, Jagst...
Japan: Ivoclar Vivadent K.K., 1-28-24-4F Hongo Bunkyo-ku Tokyo 113-0033...
Mexico: Ivoclar Vivadent S.A. de C.V., Av. Insurgentes Sur No. 863, 03810 México, D.F...
Poland: Ivoclar Vivadent Sp. z o.o., ul. Jana Pawła II 78, 01175 Warszawa...
USA: Ivoclar Vivadent Inc., 175 Pineview Drive, Amherst, NY 14228

ISO 13485 Quality Management System Certified Made in U.S.A. www.ivoclarvivadent.com Rev. Date 1/17 Rx ONLY For dental use only

INSTRUCTIONS FOR USE

MODELLATION Design the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.
SPRUNG Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.
INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.
BURN-OUT The suggested burnout temperature: 750-820C/1380-1510F
MELTING AND CASTING Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Use and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.
Casting Temperature: 1335-1395C/2435-2545F
METAL PREPARATION Carefully divest and clean the object with Al2O3. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!
OXIDATION Blast surface with 50-100 micron Al2O3 at max. 5.2 bar/75 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 1010C/1850F without vacuum with 10 min. hold time at top temperature.
HEAT TREATMENT Hardening: 540C/1000F for 15 minutes; air cool.
SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.
Pre-Solder: Special High Fusing White Ceramic Solder
Post-Solder: .615 Fine Solder, Low Fusing White Gold Solder, .585 Fine Solder, Bondal Flux
Laser weld material: Laser Ceramic Yellow
POLISHING After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.
INDICATIONS Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.
CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.
SIDE EFFECTS In individual cases, sensitivity or allergies to components of this alloy may occur.
INTERACTIONS Galvanic effects may occur between different alloys in the same oral environment.
For additional information look into the alloy property chart.

SPRUNG METHOD

RECOMMENDATIONS DIRECT: single crowns, inlays and onlays INDIREKT: multiple units and multiple single crowns
INSTRUCTIONS: 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration. 2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides. 3. Make sure that the sprue is connected to the thickest part of the restoration. 4. The connection between the sprue and the restoration should be flared ("trumpet" configuration), eliminating turbulence of the alloy (causing investment erosion) and facilitating an undisturbed flow of alloy during casting and solidification. 5. Maintain proper alloy weight in order to eliminate the negative effect of a too large button during solidification. The rule of thumb for proper alloy weight determination is: total wax weight x specific gravity of the alloy = proper alloy weight.
SUGGESTIONS: 1. Thickness and design of wax pattern: follow your manufacturer's instructions. 2. Use chill-vents when casting heavy and/or large restorations.

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