## INSTRUCTIONS FOR USE

**MODELLATION**  
Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

**SPRUNG**  
Provide the modelled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

**INVESTMENT**  
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.

**BURN-OUT**  
The suggested burnout temperature: 750-820C/1380-1510F

**MELTING AND CASTING**  
Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

**Casting Temperature:** 1285-1345C/2345-2455F

**METAL PREPARATION**  
Carefully divest and clean the object with Al<sub>2</sub>O<sub>3</sub>. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

**OXIDATION**  
Blast surface with 50-100 micron Al<sub>2</sub>O<sub>3</sub> at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 950C/1740F without vacuum with 5 min. hold time at top temperature.

**HEAT TREATMENT**  
**Hardening:** 600C/110F for 30 minutes; air cool.  
**SOLDERS AND FLUXES**  
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

**Pre-Solder:** Special High Fusing White Ceramic Solder, High-Fusing Bondal Flux  
**Post-Solder:** Low Fusing White Gold Solder, Bondal Flux  
**Laser weld material:** Laser Ceramic White

**POLISHING**  
After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

**INDICATIONS**  
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**CONTRAINDICATIONS**  
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

**SIDE EFFECTS**  
In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
Galvanic effects may occur between different alloys in the same oral environment.  
**For additional information look into the alloy property chart.**

**INDICATIONS**  
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**CONTRAINDICATIONS**  
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

**SIDE EFFECTS**  
In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
Galvanic effects may occur between different alloys in the same oral environment.  
**For additional information look into the alloy property chart.**

**INDICATIONS**  
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**CONTRAINDICATIONS**  
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

**SIDE EFFECTS**  
In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
Galvanic effects may occur between different alloys in the same oral environment.  
**For additional information look into the alloy property chart.**

**INDICATIONS**  
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns, implant superstructures, partial dentures.

**CONTRAINDICATIONS**  
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

**SIDE EFFECTS**  
In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
Galvanic effects may occur between different alloys in the same oral environment.  
**For additional information look into the alloy property chart.**

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