Instructions for Use

■ PRODUCT DESCRIPTION

Pd-based dental metal-ceramic alloy, Type 4

■ INDICATIONS*

Crowns, Telescope Crowns, Conus Crowns, Bridges, Wide Bridges, Cast Posts-/Cores, Bars, Attachments, Implant-retained Superstructures, Partial Dentures

■ WAXING/MODELLATION

Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a minimum thickness of 0.3 mm. Abutment crowns require a minimum thickness of 0.5 mm. All other cast crowns and bridges require a minimum thickness of 0.7 mm. The design of the framework provides adequate support for the veneering material. Avoid sharp angles. Connectors must have the required dimensions to provide resistance to deformation. Create large surface areas for planned soldering, with a gap of 0.05–0.2 mm.

■ SPRING

Provide the modeled single-tooth restoration or bridge framework with springs of a suitable size. In general the reservoir, sprue leads, and connector spurs, whether pear shaped or traditional, must be sized according to the specific technique. When using the direct or indirect technique be sure that the reservoir is positioned directly behind the connector between the reservoir and the casting should be a maximum of approximately 2.5–3.0 mm in length and width. The wax pattern including the spurs must be weighed in grams in order to determine the needed amount of alloy. Wax conversion formula: wax weight (gram) x alloy density = grams of alloy required

■ INVESTING

Use a phosphate-bonded investment material. Follow the manufacturer's instructions.

■ PREHEATING / BURN-OUT

Recommended burn-out temperature: 900 °C

■ MELTING AND CASTING

Flame: Propano 0.35 bar, Oxygeno 0.7 bar
Other specifics may be required by the type of casting machine. It is recommended to use a separate and clean ceramic crucible for each alloy. Preheat the ceramic crucible in the burnout furnace. The recommended ratio of used material to new material is 1:1. Do not use flux.

■ FRAMEWORK FINISHING

After being cooled, carefully dress and clean the casting with aluminum oxide (Al₂O₃). Do not use a hammer for dressing. Finishing the casting with carbide burs and/or with ceramic-bonded grinding instruments. Abrast the surface with 50–110 micron aluminum oxide (Al₂O₃) at 2.0 bar. Subsequently, steam clean or ultrasonic clean with distilled water or ethanol and dry the framework.

■ OXIDATION

Place the framework on the firing tray providing adequate support. To achieve a uniform result follow the oxidation cycle.
Temperature: 930 °C. **Holding time:** 5 min; **Vacuum:** No
If the oxide layer is stained, grind and blast the surface again. Repeat the oxide firing. Use the appropriate ceramic veneering material, following the manufacturer's instructions.
Use a ceramic Metabond® veneering material when recommended in the literature.
Höchste empfohlene Brenntemperatur: 1000 °C

■ HEAT TREATMENT

Hardening: 600 °C for 15 min; bench cool

■ SOLDERING AND LASER WELDING

The soldering gap should not be wider than the thickness of the soldering material. Follow the soldered casting to the manufacturer's instructions.
Pre Solder: Porta IP Lot-V1 **Flux:** High Fusing Bondal Flux
Post Solder: Porta OP Lot-W-2 **Flux:** Bondal Flux
Laser Welding Wire: Porta P6

■ POLISHING

Carefully remove any oxide and flux residue. Smooth the metal surfaces with rubber polishers. Polish to a high gloss finish using polishing paste. Subsequently, clean with ultrasonic cleaning equipment or careful steam cleaning.

ADDITIONAL SAFETY CONCERNS AND INSTRUCTIONS

■ CONTRAINDICATIONS

Patient with known allergy/sensitivity to any major or minor elements of this alloy, consultation with a physician is recommended. Alloy is not to be used for any application not included within the indications.

■ SIDE EFFECTS

In individual cases, sensitivity or allergies to elements of this alloy may occur. Ivoclar Vivadent makes no claims regarding the MRI-compatibility of its dental alloys. It is recommended that the patient be made aware of the potential for dental alloys to affect MRI results and to disclose the presence of dental alloys to the MRI Technician prior to conducting a test.

■ INTERACTIONS

Galvanic effects may occur between different or dissimilar alloys in the same oral environment.

■ CAUTION

Metal vapors and metal dust are harmful if inhaled. Therefore, the use of extraction equipment and / or suitable protective masks is advised!

■ STORAGE CONDITION

Store in a dry environment at room temperature.

■ DISCLAIMER

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products and their suitability and use for any purpose not explicitly stated in the Instructions. These regulations also apply if the materials are used in conjunction with products of other manufacturers.

PROCESSING DATA		
Investment material:	phosphate-bonded	
Preheating/Burn-out Temperature:	900 °C	
Crucible:	Ceramic crucible	
Casting Temperature:	1340 °C	
Oxidation:	Temperature: 930 °C; Holding time: 5 min; Vacuum: No	
CTE:	(25–500 °C): 14.2 x 10 ⁻⁶ /K	
Recommended Ceramic Material:	IP5 Style [®] , IPS InLine [®]	
Hardening:	600 °C for 15 min; bench cool	
Pre Solder / Flux:	Porta IP Lot-V1 High Fusing Bondal Flux	
Post Solder / Flux:	Porta OP Lot-W-2 Bondal Flux	
Laser Welding Wire:	Porta P6	

TECHNISCHE DATEN (ISO 22674:2016 & ISO 9693-1:2012)		
Type / Color:	4 White	4
Density (g/cm ³):	10.7	10.7
Melting Range (Solidus/Liquidus):	1100–1190 °C	1100–1190 °C
Elastic Modulus (GPa):	122	122
	Porcelain Fired	Hardened
Vickers Hardness:	410 450	410 450
Tensile Strength (MPa):	1060 1100	1060 1100
0.2% Proof Stress (MPa):	910 990	910 990
Elongation (%):	4 2	4 2

* See TYPE CLASSIFICATION DUE TO PHYSICAL PROPERTIES

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