Multilink Speed

English

Instructions for Use

Mode d'emploi

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Instructions for Use

Description

Multilink Speed is a self-adhesive, self-curing dental resin cement with light-curing option for the cementation of indirect restorations made of high-strength all-ceramic, composite, metal-ceramic and metal. When Multilink Speed is used, additional dentin bonding agents are not required. The automix syringe enables a time-saving direct application of the luting material into the restoration.

Shade

Multilink Speed is available in the shade Translucent.

Working and setting times

The following times apply as soon as Multilink Speed has been dispensed from the automix syringe (The working and setting times depend on the ambient temperature):

	At room temperature 23 °C ± 1 °C	Intra-orally 37 °C ± 1 °C
Working time	150 ± 30 sec	110 ± 30 sec
Setting time (without working time)	240 ± 60 sec	160 ± 40 sec

Attention

After Multilink Speed has been dispensed from the automix syringe, the setting process can be accelerated by intensive illumination (operating light, ambient light).

Composition

The monomer matrix is composed of dimethacrylates and acidic monomers. The inorganic fillers are barium glass, ytterbium trifluoride, co-polymer and highly dispersed silicon dioxide. Additional contents are catalysts, stabilizers and colour pigments (< 1%). The primary particle size of the inorganic fillers is between 0.1µm and 7µm. The mean particle size is 5µm. The total content of inorganic fillers is approx. 40 vol%.

Indication

- Permanent cementation of indirect restorations made of
- metal and metal-ceramic (inlays, onlays, crowns, bridges, endodontic posts)
 - high-strength all-ceramic: zirconium oxide, lithium disilicate and aluminium oxide ceramics (crowns, bridges, endodontic posts)
 - fibre-reinforced composite (endodontic posts)

6330371209WEG

CE 0123

ivoclar
vivadent®
clinical

Contraindication

Multilink Speed is contraindicated

- in situations where the preparation does not provide adequate retention (e.g. veneers).
- if a patient is known to be allergic to any of the ingredients of Multilink Speed.
- In general, Multilink Speed should not be applied to exposed pulp or to dentin which is close to the pulp.

Side effects

Systemic side effects are not known to date. In rare cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, the application of materials containing these substances is to be avoided. Only use eugenol-free temporary cements (e.g. System®.link), as eugenol may inhibit the polymerization of Multilink Speed.

Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Thus the tooth substance must not be disinfected using oxidative agents. Alkaline jet media (e.g. Airflow) impair the effect of Multilink Speed.

Application

For details please refer to the individual Instructions for Use of the products mentioned.

1. Removal of the temporary restoration and thorough cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or preparation with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray.

Then lightly dry with water- and oil-free air. Avoid overdrying.

Note:

Cleaning with alcohol can lead to overdrying of the dentin.

2. Try-in of the restoration

Next, the shade, the accuracy of fit and the occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture. If necessary, make corrections with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces. It is essential to ensure reliable isolation of the operative area with e.g. OpraGate, cotton rolls or saliva ejector. Dental hard tissue contaminated by saliva needs to be cleaned again as indicated in section 1.

3. Cleaning and pre-treatment of the restoration

3.1 High-strength glass-ceramics (e.g. IPS e.max® Press, IPS e.max CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic Etching Gel) for 20 seconds (IPS e.max Press, IPS e.max CAD). For all other products, observe the instructions of the manufacturer!
- Thoroughly rinse the restoration with water spray and dry with oil-free air.

- Apply Monobond-S or any other silanizing agent to the pre-treated surfaces with a brush or microbrush and allow it to react for 60 seconds. Subsequently, disperse Monobond Plus with a strong stream of air.
- 3.2 Clean and pre-treat the restorations made of other materials according to the instructions of the manufacturer.

Note:

- **IMPORTANT!** In order to achieve a reliable bond, do not clean zirconium oxide and base metal surfaces with phosphoric acid.
- Contamination with saliva, blood or water must absolutely be avoided during the entire procedure.
- In case of contamination, the entire pre-treatment of the restoration is to be repeated (see section 3).

4. Application of Multilink Speed into the restoration

For each application, place a new automix tip on the double syringe. Extrude Multilink Speed from the automix syringe and apply the desired quantity directly into the restoration. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until needed once again (replace with a new tip just before the next use).

5. Placement of the restoration and removal of excess cement

a) solely self-curing

Seat the restoration and retain it in place exerting uniform pressure.

Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).

b) self-curing with additional light-activation to tack-cure excess material or accelerate the curing process (quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridges)

Seat the restoration and fix/hold it in place exerting uniform pressure. Excess material is light-cured with the polymerization light (approx. 650 mW/cm², e.g. bluephase®, LOP mode) for **1 second per quarter surface** (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0-10mm. Subsequently, it can be easily removed with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). After complete curing, the removal of excess can be difficult in these areas. Subsequently, light-cure all cement joints again for 20 seconds (approx. 800 mW/cm²; e.g. bluephase C8, HIP mode) for final polymerization.

Note

As with all composites, Multilink Speed is subject to oxygen inhibition. This means that the surface layer does not polymerize during curing, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

6. Finishing of the completed restoration

- Apply finishing and polishing strips in the proximal regions.
- Check occlusion and functions and adjust if necessary.
- Polish restoration margins with silicone polishers (OpraPol, Astropol®, Politip®-F) or discs.

Special note

Multilink Speed should have room temperature when it is applied. Cool temperatures can render the material difficult to extrude and increase the processing and curing times.

Warning

Unpolymerized Multilink Speed paste is an irritant. Avoid contact with the skin, mucous membrane and eyes. If Multilink Speed comes into contact with the eyes, immediately rinse with copious amounts of water and seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use Multilink Speed after the indicated date of expiration.
- Cool storage (2–8 °C/36–46 °F) for Multilink Speed.
- Store Multilink Speed automix syringe with mixing tip attached after use.
- Shelf life: see expiration date.

Keep material out of the reach of children!

For use in dentistry only!

Date information prepared: 12/2009, Rev 0

Manufacturer:

Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

The Material Safety Data Sheet is available online at

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.