**Descriptions**

- **N-Etch** is a phosphoric acid gel for enamel etching and dentin conditioning.
- **Tetric® N-Bond** is a light-curing, single-component bonding agent for enamel and dentin bonding in combination with the total-etch technique.
- **Tetric® N-Flow** is a flowable, light-curing, radiopaque nano-hybrid composite for restorative treatment and the cementation of translucent ceramic and composite restorations.
- **Tetric® N-Ceram** is a light-curing, radiopaque nano-hybrid composite for restorative treatment. Tetric N-Ceram, Tetric N-Flow and Tetric N-Bond cure with light in the wavelength range of 400–500 nm.

**Shades**

**Tetric N-Ceram** is available in the following shades:

*Enamel shades:* A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D3

*Dentin shades:* A2 Dentin, A3.5 Dentin, B2 Dentin

*Highly translucent Incisal shade:* T

*Bleach shades:* Bleach Light (L), Bleach Incisal (I), Bleach Medium (M)

**Tetric N-Flow** is available in the following shades:

*Enamel shades:* A1, A2, A3, A3.5, A4, B2

*Dentin shades:* A3.5 Dentin

*Highly translucent Incisal shade:* T

*Bleach shades:* Bleach Light (L), Bleach Incisal (I)

**Composition**

- **N-Etch** contains phosphoric acid (37 wt.% in water), thickeners and pigments.
- **Tetric N-Bond** contains phosphoric acid acrylate, HEMA, Bis-GMA, urethane dimethacrylate, ethanol, film-forming agent, initiators and stabilizers.
- **Tetric N-Flow** contains 36 wt.% dimethacrylates (including TEGDMA), 63 wt.% fillers (barium glass, ytterbium trifluoride, highly dispersed silica and mixed oxide) and 1 wt.% initiators, stabilizers and pigments. The total content of inorganic fillers is 39 vol.%. The particle size of inorganic fillers is between 40 and 3000 nm.
- **Tetric N-Ceram** consists of dimethacrylates (19-20 wt.%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (80-81 wt.%). Additives, initiators, stabilizers and pigments are additional contents (< 1 wt.%). The total content of inorganic fillers is 55–57 vol.%. The particle size of inorganic fillers is between 40 nm and 3000 nm.

**Indications**

- **N-Etch** is a gel used for enamel etching and dentin conditioning in conjunction with adhesives, composites, luting materials and fissure sealants.
- **Tetric N-Bond**
  - Adhesive for direct, light-curing and dual-curing composite and compomer restorations
  - Adhesive for indirect all-ceramic and composite restorations (inlays, onlays, veneers)
- **Tetric N-Flow**
  - As initial layer under Class I and II restorations
  - Anterior restorations (Classes III, IV)
  - Class V restorations (cervical caries, root erosion, wedge-shaped defects)
  - Small restorations of all types
  - Extended fissure sealing
  - Splinting of mobile teeth
  - Blocking out of undercuts
  - Adhesive cementation of translucent ceramic and composite restorations
- **Tetric N-Ceram**
  - Class I-V restorations
  - Restoration of deciduous teeth
  - Extended fissure sealing
  - Direct veneers
  - Splinting of mobile teeth
  - Blocking out of undercuts
  - Repair of composite and ceramic veneers

**Contraindications**

- **N-Etch**
  - If the stipulated application technique cannot be applied.
- **Tetric N-Bond**
  - If the patient is allergic to any of the components of Tetric N-Bond or if the stipulated working technique cannot be used.
  - Tetric N-Bond must not be used in combination with purely self-curing composite materials. Always activate dual-curing materials with light.
  - Direct pulp cappings
- **Tetric N-Flow/Tetric N-Ceram**
  - If a dry working field cannot be established or the stipulated application technique cannot be applied.
  - If the patient is known to be allergic to any of the ingredients.

**Side effects**

- In rare cases, components of Tetric N-Bond/Tetric N-Flow/Tetric N-Ceram may lead to sensitization. In such cases, the products should no longer be used.
- To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp).

**Interactions**

- Glass ionomer-based liners may be dissolved by N-Etch.
- Substances containing eugenol/clove oil may inhibit the polymerization of adhesives and composite materials. Consequently, the application of such materials together with Tetric N-Bond/Tetric N-Ceram/Tetric N-Flow must be avoided.
- Cationic mouth rinses and plaque disclosing agents as well as chlorhexidine may cause discolouration.
- Materials containing peroxides (e.g. peroxide-based bleaching agents or disinfectants) may inhibit the polymerization of adhesives.
- Phenolic substances (e.g. eugenol, salicylic derivatives) inhibit polymerization of Tetric N-Ceram/Tetric N-Flow/Tetric N-Bond.
- Iron- und aluminium-based styptic agents must not be used shortly before or during the application of Tetric N-Bond.
- Fluoridation agents must not be used prior to the application of Tetric N-Bond.

# Tetric® N-Ceram Tetric® N-Flow Tetric® N-Bond N-Etch

**EN Instructions for Use**  
Light-curing dental restorative composite, flowable composite, dental adhesive and etching gel

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## Application

### Direct Restorations

- Select the shade  
For facilitated access to the treatment field, the lip and cheek retractor OptraGate® can be used.
- Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.
- Prepare the cavity according to the requirements of the adhesive technique.
- Clean the cavity with water spray.
- Air dry the cavity.
- Apply a pulp protector (calcium hydroxide material, e. g. ApexCal®) if needed; only cover areas close to the pulp and subsequently apply a pressure resistant cement (e.g. Vivaglass® Liner).
- Place a matrix and an interdental wedge, if required.
- Apply N-Etch on enamel and subsequently on dentin; and allow a reaction time of 15 seconds. Afterwards, thoroughly rinse off the etchant with water spray and dry the tooth surfaces with oil-free air. Avoid excessive drying of the dentin. The etched enamel surface should have a chalky white appearance. If this is not the case, repeat the etching procedure using a shorter reaction time. If the tooth surface becomes accidentally contaminated (e. g. with saliva), repeat the etching procedure with a reaction time of max. 10 seconds. The reaction time on unprepared enamel (e.g. fissure sealants) is 30-60 seconds. If the enamel is selectively etched, use a reaction time of 15-30 seconds.
- Tetric N-Bond.  
Dispense the desired amount of Tetric N-Bond into a mixing well and apply it using a disposable applicator (e.g. Vivadent® Applicator Brush). Protect the adhesive from light (e.g. VivaPad®). Close the bottle immediately after each use. A new disposable applicator must be used for each application. Apply a thick layer of Tetric N-Bond on the enamel and dentin surfaces, using the enclosed application brush. Brush the material gently into the dentin for at least 10 seconds. Avoid insufficient coverage of the cavity and apply additional material, if required. Remove excess material and the solvent by a gentle stream of air so that the adhesive completely covers the enamel and dentin without pooling.
- Light-cure Tetric N-Bond according to the curing time recommendations below.  
A shiny tooth surface prior to the application of the composite shows that all surfaces are completely covered.
- Apply Tetric N-Flow in layers of max. 2 mm or 1.5 mm (Dentin shade) and light-cure according to the curing recommendation indicated below. Apply Tetric N-Ceram in layers of max. 2 mm or 1.5 mm (Dentin shade) and contour/adapt the material to the cavity walls using a suitable instrument (e.g. OptraSculpt®).

Polymerize each layer individually according to the table below.

Hold the light emission window as closely as possible to the surface of the restorative material.

Light intensity	Tetric N-Ceram	Tetric N-Flow
≥ 500 mW/cm <sup>2</sup>	20 s	20 s
≥ 1000 mW/cm <sup>2</sup>	10 s	10 s

Curing Light	Material	High Power program	Soft Power program	Low Power program
Bluephase N	Tetric N-Bond	–	–	10 s
	Tetric N-Flow	10 s	15 s	–
	Tetric N-Ceram	10 s	15 s	–
Bluephase NM	Tetric N-Bond	10 s	–	–
	Tetric N-Flow	15 s	–	–
	Tetric N-Ceram	15 s	–	–
Bluephase NMC	Tetric N-Bond	10 s	–	–
	Tetric N-Flow	15 s	–	–
	Tetric N-Ceram	15 s	–	–

- Finish the restoration with suitable finishers or fine diamonds. Subsequently, check the occlusion. Polish with silicone polishers (e.g. OptraPol®).

### Indirect Restorations

- Prepare the ceramic or composite surface according to the manufacturer's instructions.
- Prepare the tooth in the same way as for direct restorations, steps 1-10.  
**Important notes:** Avoid pooling of the adhesive as this may compromise the accuracy of fit of the definitive restoration.  
Light-cure the adhesive before placing the restoration.

- Subsequently, place the restoration using a light-/dual-curing luting material (e.g. Variolink N) and light-cure according to the instructions of the manufacturer.

### Note on N-Etch Jumbo

#### Refilling of empty N-Etch syringes

Empty 2 ml syringes can be conveniently refilled with the contents of the N-Etch Jumbo syringe up to 15 times.

#### Important

- Although the 2 ml syringe is equipped with a plunger safety stop, it is absolutely necessary to wear gloves, protective goggles and a face mask for reasons of safety.
- Remove the screw-on cap from the 30 g Jumbo syringe and the empty 2 ml syringe.
- Join the tips of the two syringes with a twisting, clockwise motion. Make sure that the two syringes are firmly joined.  
**Caution:** Do not strip the thread.
- Refill the 2 ml syringe by slowly and steadily pulling the plunger of the small syringe. Simultaneously, exert light pressure on the plunger of the jumbo syringe using the other hand to support the refilling process.
- Avoid overfilling the 2 ml syringe.**
- Under no circumstances exceed the plunger safety stop of the small syringe – there is a risk of acid burn!**
- Do not point the syringes towards people during refilling.
- Do not fill etchant from the 2 ml syringe back to the Jumbo syringe.
- Disconnect the two syringes by an anticlockwise twist.
- Replace the screw-on cap on the 30 g Jumbo syringe.
- Attach an application tip to the 2 ml syringe. The filled 2 ml syringe is now ready for use. Apply the etchant as described above.
- Use only the N-Etch application tips delivered with this package. The use of other application tips may have an adverse effect on the viscosity of the etchant.
- In order to ensure traceability, write the LOT number and expiration date of the 30 g Jumbo syringe on the label provided; and attach it to the filled 2 ml syringe.

#### Special note

- If fresh Tetric N-Ceram is applied to polymerized material, it must first be roughened and wetted with Heliobond.
- Tetric N-Flow/Tetric N-Ceram should have room temperature when applied. Cool temperatures render the material difficult to extrude.
- Do not use oxidizing disinfectants to disinfect bottles, syringes and Cavifils.

#### Warning N-Etch

Phosphoric acid gel is corrosive. When filling the 2 ml syringe, always wear protective goggles and gloves. Avoid contact with the skin, mucous membrane and eyes (goggles are recommended for both patients and operators). If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician/ophthalmologist. After skin contact, immediately rinse with water and soap. After contact with sensitive materials (surfaces, clothing), immediately clean with copious amounts of water.

#### Warning Tetric N-Bond

- Tetric N-Bond is an irritant. Avoid contact with skin, mucous membranes and eyes.
- In case of skin contact, immediately wash with soap and copious amounts of water.
- If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician/ophthalmologist.
- In rare cases, contact with skin may lead to sensitization towards the ingredients.
- Commercial medical gloves do not provide protection against the sensitizing effects of methacrylates.

#### Warning Tetric N-Flow/Tetric N-Ceram

- Avoid contact of unpolymerized material with the skin/mucous membrane and eyes.
- Unpolymerized material may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

#### Storage and hygiene

- Storage temperature: 2-28 °C (36-82 °F)  
N-Etch recommendation: Refrigerate large stocks at 2-8 °C (36-46 °F).
- Close syringes/Cavifils/bottles immediately after usage.  
Exposure to light causes premature polymerization.
- Expiry date: see note on Cavifils, syringes, bottles and packaging.
- Do not use any of these products after the indicated date of expiration.
-  If Tetric N-Ceram or Tetric N-Flow is applied from the Cavifil directly in the mouth of the patient, the Cavifil or the application cannula of the syringe must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

#### Keep material out of children's reach.

#### For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Intensité lumineuse	Tetric N-Ceram	Tetric N-Flow
≥ 500 mW/cm <sup>2</sup>	20 sec.	20 sec.
≥ 1000 mW/cm <sup>2</sup>	10 sec.	10 sec.

Lampe à photopolymériser	Matériau	Programme High Power	Programme Soft Power	Programme Low Power
Bluephase N	Tetric N-Bond	-	-	10 sec.
	Tetric N-Flow	10 sec.	15 sec.	-
	Tetric N-Ceram	10 sec.	15 sec.	-
Bluephase N M	Tetric N-Bond	10 sec.	-	-
	Tetric N-Flow	15 sec.	-	-
	Tetric N-Ceram	15 sec.	-	-
Bluephase N MC	Tetric N-Bond	10 sec.	-	-
	Tetric N-Flow	15 sec.	-	-
	Tetric N-Ceram	15 sec.	-	-

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#### Tetric N-Flow/Tetric N-Ceram - Avertissement

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**Ivoclar Vivadent Marketing Ltd.**
**Ivoclar Vivadent AB**
**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**