Description

Systemp®.c&b II is a self-curing composite material used for the fabrication of temporary inlays, onlays, crowns and bridges. It is supplied in specially designed cartridges that provide a 10:1 base/catalyst ratio.

Composition

Polyfunctional methacrylates (48 wt.%), inorganic fillers (47 wt.%), plasticizers, initiators, stabilizers and pigments (5 wt.%).

Indications

Systemp.c&b II is a temporary crown and bridge material suitable for the fabrication of inlays, onlays, crowns and bridges. In addition, Systemp.c&b II may be used for lining prefabricated polycarbonate crowns.

Contraindication

Do not use if a patient is known to be allergic or intolerant to any of the ingredients of Systemp.c&b II.

Side effects

In rare cases, components of Systemp.c&b II may lead to a sensitizing reaction. In this case, the use of the product should be discontinued.

Application**1. Fabrication of preliminary impression****(Alginate or silicone matrix):**

If the unprepared teeth exhibit a suitable occlusal morphology that may be transferred to the provisional restoration or even the planned permanent restoration, and if abutment teeth in suitable position are available, a preliminary impression can be made to serve as a matrix for the fabrication of the temporaries. Prior to the preparation or extraction of the teeth, an impression of the preoperative situation is made using addition cross-linked silicone impression material (e.g. Virtual®; dimensionally stable impressions) or alginate (impressions that are dimensionally stable for a limited period of time). To facilitate repositioning of the impression in the oral cavity, the interdental gingival areas should be trimmed and undercuts removed in the impression, if necessary. If the occlusion is defective, the fabrication of a wax-up in the laboratory may additionally be required prior to taking the impression.

Plastic template:

When using a vacuum-formed plastic matrix, care should be taken that the plastic material does not bond chemically to Systemp.c&b II. Suitable plastic sheets (e.g. polyethylene sheets) should be selected, so that the vacuum-formed matrix can be reused. If a bond between the plastic and Systemp.c&b II is desired, plastic sheets that are especially designed for this purpose should be chosen. The interdental aspect of the matrix can be conditioned with a bonding agent (e.g. Heliobond) prior to applying the temporary restoration composite for adhesion.

2. Preparing the double cartridge

Fig. 1

1. Press the black release lever located below the plunger on the back of the dispenser and pull the plunger as far back as possible (Fig. 1).

Fig. 2

2. Lift the cartridge lock and insert the cartridge. Lower the cartridge lock (Fig. 2). Push slide forward until it touches the cartridge.

Fig. 3

3. Remove the sealing cap or used mixing tip by turning it a 1/4 turn counter clockwise and discard it, do not reuse it (Fig. 3).

Fig. 4

4. Insert a new mixing tip. Push the tip down until the notch on the mixing tip touches the notch on the cartridge. While gripping the coloured base, and not the mixing tip, turn the mixing tip base 1/4 turn clockwise (Fig. 4).

5. The material is mixed and dispensed through the mixing tip by pulling the trigger and can be applied directly into the matrix.

Removal of the cartridge

Lift the release lever below the plunger on the back of the dispenser and pull the plunger as far back as possible. Lift the cartridge lock on the front side on top of the dispenser and remove the cartridge.

Leave the used mixing tip on the syringe as a cap until the next application.

3. Application

Systemp.c&b II is automatically mixed when it is dispensed with slight pressure directly into the impression or vacuum-formed matrix.

Prior to each application, extrude and discard a small amount (approx. 100 mg) of material, see Fig. 5.

To prevent voids, apply Systemp.c&b II to the occlusal surfaces first, then proceed by filling the gingival areas.

Keep the mixing tip immersed in the material to prevent bubbles. If required, Systemp.c&b II may also be syringed

around the prepared teeth to avoid bubble formation along the preparation margin. The processing time is approx. one minute at room temperature (23°C / 73°F).

4. Shaping of temporaries**a) Fabrication in the mouth**

The prepared teeth should be slightly moist. This can be achieved by painting them with water-soluble glycerine gel (e.g. Liquid Strip).

Undercut areas of the tooth preparation or adjacent teeth may have

Systemp.c&b II

EN	Instructions for Use Self-curing temporary c&b material	PT	Instruções de Uso Material autopolimerizável para coroas e pontes provisórias
FR	Mode d'emploi Matériau autopolymérisant pour couronnes et bridges provisoires	TR	Kullanma Talimatı Kendiliğinden sertleşen geçici kârun ve köprü materyali
ES	Instrucciones de uso Material provisional C&B autopolimerizable	RU	Инструкция по применению Самотвердеющий композитный материал для временных коронок и мостов

to be blocked out with appropriate materials (e.g. with wax). Carefully reposition the impression filled with Systemp.c&b II on the prepared teeth. The setting time in the mouth (37 °C / 98 °F) is approximately 1 to 2 minutes. By then, Systemp.c&b II has attained a hard-elastic consistency and can be removed from the mouth together with the impression.

b) Fabrication on the model

Isolate the preparation area of the model (e.g. using Vaseline). Carefully reposition the impression filled with Systemp.c&b II on the model. After approx. 3 minutes (at 23 °C / 73 °F), Systemp.c&b II attains a hard-elastic, partially-set consistency and can be removed from the model together with the impression.

5. Curing and finishing

After removing the temporary restoration from the impression (or from the tooth preparation), remove excess material using rotary instruments once the material is cured completely (after approx. 4 to 5 minutes). Suitable finishing instruments are cross-cut tungsten carbide burs. Remove the oxygen-inhibited layer with e.g. alcohol or by polishing the surface of the restoration. For polishing use silicon carbide rubber polishers (e.g. Astropol®).

6. Cementation of temporaries

In order to reduce the sensitivity of the prepared dentin surfaces during the temporization phase, Systemp®.desensitizer may be applied prior to cementing the temporary restoration. For this purpose, dentin surfaces must be clean and dry. Brush Systemp.desensitizer into the dentin surfaces for 10 seconds using a suitable instrument (brush, application aid). Cautiously disperse excess to a thin layer with a stream of air. Subsequently seat the temporary restoration, preferably using a non-eugenol temporary cement (e.g. Systemp®.link).

7. Relining, repairs, add-ons

If Systemp.c&b II temporaries are relined, repaired or adjusted, the following procedure is recommended:

Systemp.c&b II temporaries can be repaired with freshly mixed Systemp.c&b II after having been conditioned with Heliobond.

1. Grind the areas of the Systemp.c&b II temporary which need to be relined, repaired or adjusted with a coarse diamond.
2. Apply Heliobond.
3. Disperse Heliobond with a very weak stream of air. Avoid pooling! Care should be taken not to blow away the bonding agent. All the prepared surfaces must be adequately coated with Heliobond.
Note: Heliobond does not contain a solvent that has to be evaporated!
4. Light cure Heliobond for 10 seconds with a light intensity of 500 to 1,100 mW/cm² (e.g. Bluephase® N).
5. Apply Systemp.c&b II.
6. If necessary, finish the relined, repaired or adjusted areas of the Systemp.c&b II temporary with e.g. cross-cut tungsten carbide burs and subsequently polish them with silicon carbide rubber polishers (e.g. OpraPol®).

As an option, Systemp.c&b II temporaries can be directly repaired or adjusted using a light-curing composite (e.g. Tetric® N-Flow). In this case, the respective Instructions for Use must be observed.

Working time

Approx. 1 minute at room temperature (23°C / 73°F).

Setting time

Approx. 1 to 2 minutes at 37 °C / 98 °F, approx. 3 minutes at room temperature (23 °C / 73 °F).

Additional information

- Uncured material can be removed using alcohol or similar solvents.

Warnings

- Prevent any contact of uncured Systemp.c&b II with the skin, mucous membrane and eyes.
- Uncured Systemp.c&b II may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Shelf life: see expiration date on label
- Storage temperature: 12 – 28 °C / 54 – 82°F.

- Only use Systemp.c&b II at room temperature (cooled material is more viscous and exhibits a slower setting behaviour).
- Store cartridge in a horizontal position or in the dispenser, do not store in an upright position.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.