**Description**

Amalcap® Plus is a non-gamma-2, spherical / irregular-particle silver amalgam in self-activating capsules for dental restorations.

**Composition**

**One Size-1 Amalcap Plus non-gamma-2 capsule contains the following active substances:**

	regular	fast
<b>Powder</b>		
Silver	294 mg	302 mg
Tin	76 mg	77 mg
Copper	50 mg	51 mg
<b>Liquid</b>		
Mercury	400 mg	400 mg

**One Size-2 Amalcap Plus non-gamma-2 capsule contains the following active substances:**

	regular	fast
<b>Powder</b>		
Silver	424 mg	434 mg
Tin	109 mg	112 mg
Copper	72 mg	74 mg
<b>Liquid</b>		
Mercury	590 mg	590 mg

**One Size-3 Amalcap Plus non-gamma-2 capsule contains the following active substances:**

	regular	fast
<b>Powder</b>		
Silver	555 mg	570 mg
Tin	142 mg	146 mg
Copper	94 mg	97 mg
<b>Liquid</b>		
Mercury	800 mg	800 mg

Length of the capsule: 31.6 mm  
Largest outer diameter: 13.6 mm

**Indications**

Class I and II restorations

**Contraindication**

Amalcap Plus non-gamma-2 is contraindicated

- if proximal or occlusal contacts with other metal restorations are present
- if the patient suffers from impaired renal function
- if the patient is known to be allergic to any of the components of amalgam
- for retrograde or endodontic restorations
- as a core build-up material under crowns or inlays
- as a filling material under metal crowns
- for children under fifteen years of age
- for pregnant or nursing women
- if the stipulated working technique cannot be applied.

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**EN Instructions for Use**

- Non-gamma-2, 70 % silver amalgam alloy with spherical / irregular particles in self-activating capsules

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## Side effects

- In individual cases, hypersensitivity, allergic reactions, or electrochemically-induced local reactions may occur.
- Mercury may have a sensitizing effect on the skin and lungs, or act as a nephrotoxin or neurotoxin.
- Electrochemical processes may cause lichen planus of the oral mucous membrane.
- After the placement or removal of amalgam restorations, a short-term increase of the mercury concentration in the blood or urine may occur.

## Interactions

Mercury that has escaped during condensing or amalgam that has not set may cause amalgamation or galvanic reactions (electrochemically-induced local malaise, i.e. alterations in the sense of taste) when the restoration comes in contact with other metal restoratives. If these symptoms persist, the amalgam should be replaced with another material.

## Application

Relative isolation with e.g. cotton rolls.

Moisture contamination of the amalgam before it has set completely can detrimentally affect physical properties such as stability and corrosion resistance.

## Mixing times for Amalcap Plus non-gamma-2

The mixing time in the Silamat® S6 is 10 s.

Mixers of other competitors should demonstrate a high starting torque and achieve approx. 4,500 rpm. The use of low-speed amalgamators is not recommended. Longer mixing times result in a softer consistency and slightly improved stability, without significant effect on the other properties. The triturated amalgam should be bright and have a plastic consistency.

## Base/liner

Cover very deep areas of the cavity near the pulp with a small amount of calcium hydroxide cement. The base/liner must be able to withstand condensation pressure. Line the cavity with a mechanically stable cement such as Vivaglass Liner (light-curing glass ionomer cement).

## Placing the restoration

After mixing, open the capsule, apply the pliable amalgam to the cavity in portions and condense and model it. Ultrasonic condensers should not be used. The working time of Amalcap Plus Fast is approx. 3 minutes; the working time of Amalcap Plus Regular is approx. 4 minutes.

## Polishing

The restoration should not be polished until at least 24 hours after placement. Suitable tungsten-carbide finishers or other polishing instruments (Politip-F and Politip-P) can be used. Make sure to provide steady, adequate spray cooling and suction.

## Safety notes

- Mercury reacts with and embrittles particular metals and their alloys. Avoid unnecessary contact between mercury and those metals (and their alloys).
- Mercury presents a health hazard if handled incorrectly. Mercury is toxic by vapour inhalation and the effect is cumulative. Spillages of mercury should be removed immediately, including from places which are difficult to access. Use a plastic syringe to draw it up. Smaller quantities can be covered by sulphur powder and removed. Individual small droplets can be picked up by tin foil (Sn) and removed. Avoid inhalation of the vapour. During and after clean-up thoroughly ventilate the area where the spill has occurred. Vacuum cleaners should not be used.
- If moisture is introduced into the amalgam before it has set, properties such as strength and corrosion resistance can be affected adversely. Use a dry field, whenever possible.
- Waste material and all primary containers that have held mercury shall be disposed of following appropriate management practice.
- As the removal of existing amalgam restorations results in the additional release of mercury, clinically sound amalgam restorations should not be removed. The removal of existing amalgam restorations should be avoided during pregnancy or lactation in particular.
- When removing amalgam restorations, the corresponding safety notes must be observed: Wear nose/mouth protection and goggles; use spray aspirators, saliva ejectors and provide adequate spray cooling.
- Only store amalgam and place and remove amalgam fillings in well-aired practice facilities.
- Avoid skin contact with unset material.
- Never open capsules prior to mixing!
- Single-use only

## Storage

- **Store amalgam waste in containers that can be tightly sealed.**
- **Dispose of used capsules in accordance with national regulations.**
- Store packages at room temperature (2–28°C/36–82°F).
- See expiration date for shelf life.

## Keep material out of children's reach!

### For dental use only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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**Ivoclar Vivadent Marketing Ltd.**
**Ivoclar Vivadent AB**
**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**