Description

Systemp.link is a dual-curing (light- and self-curing) temporary composite cement to be used for the aesthetic temporary cementation of provisional crowns, bridges, inlays, onlays and veneers. Systemp.link is compatible with all common temporary composite materials (e.g. Systemp®.c&b II).

Shades

Transparent

Composition

Systemp.link is composed of bismethacrylates (56 wt. %) and fillers (43 wt. %). Catalysts, stabilizers and pigments are additional ingredients.

Indications

- Temporary cementation of provisional restorations
- The product has been designed to remain in the mouth for up to 4 weeks.

Contraindication

The product should not be used in patients with a known allergy to any of its ingredients.

Side effects

In rare cases, components of Systemp.link may lead to a sensitizing reaction in susceptible individuals. In these cases, the use of the product should be discontinued.

Interactions

The eugenol-free formula makes Systemp.link an ideal temporary cement even if adhesive luting techniques are subsequently employed. Eugenol-containing cements are contraindicated in these cases, as they may affect the polymerization of the luting composite used for the permanent cementation of the final restoration.

To prevent Systemp.link from permanently adhering to composite-based core build-up materials, the composite core should be coated with a thin layer of separator. Water-soluble glycerine gel (e.g. Liquid Strip) may be used for this purpose. If dentin adhesives or protective varnishes are used after preparation, subsequent isolation with Vaseline or glycerine gel is required.

Application

To reduce dentinal hypersensitivity and postoperative sensitivity, Systemp.desensitizer should be applied according to the IfU of Systemp.desensitizer before temporary cementation.

Application of the Systemp.desensitizer:

- 1) The dentin surfaces should be dry and clean.
- 2) Apply Systemp.desensitizer on the dentin and brush into the dentin for approximately 10 seconds using a suitable instrument (brush, application brush).
- 3) Carefully disperse excess with blown air / dry with an air syringe.

Application of the Systemp.link:

- 1) Remove the cap of the double-push syringe by turning it 1/4 turn counter clockwise (discard the cap, do not re-use it!) and replace it with a mixing tip. If the double-push syringe has been used before, replace the used mixing tip with a new one just before dispensing Systemp.link.
- 2) When attaching the mixing tip to the syringe, make sure that the guidance of the double-push syringe is aligned with that of the mixing tip. Push the tip completely down until the notch on the mixing tip is aligned with that of the double-push syringe. Secure the mixing tip in place by gripping the coloured base and turning it 1/4 turn clockwise. The double-push syringe contains predosed amounts of Systemp.link base and catalyst which are automatically mixed and dispensed when the two components are extruded. Thus the material may be directly applied to the provisional restoration.
- 3) Spread a layer of Systemp.link over the dry inner surfaces of the provisional restoration and/or the prepared, cleaned and dried tooth surfaces. The working time is approx. 2.5 to 3 minutes at a room temperature of 23 °C / 73 °F.
- 4) Subsequently, seat the provisional restoration with firm but light pressure. Excess cement can be removed after approx. 3 minutes using a scaler or other appropriate instruments. Also use a scaler and dental floss to remove any excess cement from the interproximal areas. To accelerate the polymerization process, cure the cement layer with a high performance curing light (light intensity > 1000 mW/cm²; e.g. Bluephase® N) for 10 seconds per aspect (occlusal, mesial, distal), or 20 seconds per aspect (occlusal, mesial, distal) with a standard curing light (light intensity > 500 mW/cm²). By additionally curing the crowns from the occlusal aspect, the primary stability can be increased. The light curing procedure allows excess cement to be removed immediately.

Leave the used mixing tip on the syringe.

It will serve as a seal until the next use!

Special notes

- To ensure an easy clean-up and complete removal of excess, we recommend holding back the lip or cheek surfaces during the setting process, to prevent the excess material from spreading out into a thin layer.
- During the self-curing process, the patient should not exert any pressure on the provisional restoration.

Working time

Approx. 2.5 to 3 minutes at a room temperature of 23 °C / 73 °F.

Self-curing time

4 minutes at a temperature of 37 °C / 98 °F. Exposure to light will accelerate the curing process.

Note

- Extruding material without the mixing tip in place may result in a cross-contamination of the Systemp.link components.
- Like all other composites, Systemp.link is subject to oxygen inhibition. Therefore, the behaviour of the cement when mixed on a mixing pad

Systemp.link

EN Instructions for Use

Dual-curing resin based temporary cement

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outside the mouth cannot be taken as a reference for the intraoral setting behaviour.

- While the temporary restoration is in place, the material may undergo some change in colour. However, this does not affect the functionality.

Precautions

- Prevent any contact of uncured material with the skin/mucous membrane and eyes. Unpolymerized Systemp.link may cause slight irritation and in very rare cases, sensitization may occur if the skin is exposed repeatedly.
- Commercially available medical gloves do not offer protection against the sensitizing effect of methacrylates.

Storage and shelf life

- Shelf life: see expiration date on label/packaging
- Do not use Systemp.link after the indicated date of expiration.
- Storage temperature: 2 to 8 °C / 36 to 46 °F
- Only use Systemp.link at room temperature (cooled material is more viscous and exhibits a slower setting behaviour).

Keep out of the reach of children!

For use in dentistry only!

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Tel. +82 2 536 0714, Fax +82 2 596 0155, www.ivoclarvivadent.co.kr**Ivoclar Vivadent S.A. de C.V.**
Tel. +31 23 529 3791, Fax +31 23 555 4504, www.ivoclarvivadent.com**Ivoclar Vivadent Ltd.**
Tel. +48 22 635 5496, Fax +48 22 635 5469, www.ivoclarvivadent.pl**Ivoclar Vivadent Marketing Ltd.**
Tel. +34 913 757 820, Fax +34 913 757 838, www.ivoclarvivadent.es**Ivoclar Vivadent AB**
Tel. +90 212 343 0802, Fax +90 212 343 0842, www.ivoclarvivadent.com**Ivoclar Vivadent Limited**
Tel. +44 116 284 7880, Fax +44 116 284 7881, www.ivoclarvivadent.co.uk**Ivoclar Vivadent, Inc.**