Instructions for use

Page 2

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Seite 7

Mode d'emploi

Page 13

Istruzioni d'uso

Pagina 18

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Página 24

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Página 29

Bruksanvisning

Página 35

Brugsanvisning

Página 40

Käyttöohje

Página 45

Productinformatie

Página 50

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Instructions for Use

Description

SR Ivocron is a PMMA veneering material. Depending on the indication and the monomer used, SR Ivocron can be used in conjunction with the cold, hot, or pressing technique.

Indication

Cold technique (pouring technique)

- For short- or medium-term temporaries
- Securing of ground denture teeth on the cast framework
- Repairs

Hot technique (layering technique)

- Crown and bridge veneering technique

Pressing technique (flasking technique)

- Crown and bridge veneering technique

Contraindication

Direct use in the oral cavity

Composition

SR Ivocron	Cold Liquid	Hot Liquid % wt.	Press Liquid
Methyl methacrylate	86.5	65	90
Dimethacrylate	13.0	35	10
Catalyst	0.5	–	–

- *SR Ivocron Dentin, Incisal, Cervical, and Intensive Powders*

consist of polymethyl methacrylate (> 98 % wt.), catalysts, and pigments (< 2 % wt.).

- *SR Ivocron Opaquer Powder*

consists of copolymer, aluminium oxide, barium sulphate, and titanium dioxide (> 98 % wt.), as well as catalysts and pigments (< 2 % wt.).

- *SR Ivocron Opaquer Liquid*

Methyl methacrylate ≥ 99 % wt., catalysts ≤ 1 % wt.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to PMMA materials have been reported.

Warning

SR Ivocron Opaquer Liquid, as well as the Cold, Hot, and Press Liquids contain methyl methacrylate. MMA is irritant and highly flammable (flash point: 10 °C/50 °F). Do not inhale vapours. The material irritates eyes, respiratory organs, and skin. Skin contact may lead to sensitization. For further information, please refer to the EEC safety data sheet (or MSDS).

Cold technique (pouring technique)

With the Cold Liquid, SR Ivocron may be used for short- and medium-term temporaries, securing of ground denture teeth, and repair of PMMA.

Temporaries

Designing the model, wax-up, and silicone key

Fabricate the working model from die stone. Produce the wax-up and check shape and function. After that, a silicone key is made.

Isolation

Boil out the wax and block out undercuts with pink wax. Then, immerse the model in water for 5 minutes. Apply two layers of Separating Fluid. Wait until the surface is no longer shiny before applying the second layer.

Tip

Slightly rough or undercut dies may be isolated and blocked out with the elastic SR Ivocron Separator isolating gel. Apply SR Ivocron Separator and allow to dry for approx. 5 minutes.

Incisal and Dentin

Mix the corresponding SR Ivocron incisal material in a small, clean rubber cup, cover it, and let it sit for approx.

2–3 minutes. Use the same procedure for the Dentin

material. Apply the Incisal material in the incisal area of the silicone key, while it is still flowable. As soon as the Incisal material stops flowing, the silicone key is filled up with Dentin material. Reposition the filled silicone key on the isolated model and hold it in place with a rubber band.

Polymerization

Polymerization is carried out in the pressure apparatus (e.g. Ivomat) for 15 minutes at 40–50 °C (104–122 °F) and 2–6 bar pressure (29–85 psi).

Finishing, polishing

Finish the restoration according to the standard procedure for PMMA materials, i.e. finishing with burs and rubber polishers, pre-polishing with goat's hair brushes and pumice, and polishing to a high gloss with cotton wool buffers and polishing medium (e.g. Ivoclar Vivadent Universal Polishing Paste).

Securing of ground denture teeth on cast frameworks

Setup, fabrication of the silicone key

Grind the denture tooth to fit on the working model, hold in place with wax, and check occlusion. Fabricate a silicone key.

Isolation

Boil out the wax and block out undercuts with pink wax. Then, immerse the model in water for 5 minutes. Apply two layers of Separating Fluid. Wait until the surface is no longer shiny before applying the second layer.

Wetting

Roughen the basal areas of the teeth, replace in the silicone key, and wet with monomer.

Metal bond

Use aluminium oxide to blast the metal parts that will be covered.

Tip

We recommend using SR Link for the chemical bond to metal. Blast framework with aluminium oxide and remove grinding dust. Wet it with SR Link and allow it to dry for 3 minutes.

Opaquer

Mix the desired Opaquer with the SR Ivocron Opaquer Liquid and allow a dough time of approx. 2–3 minutes. After that, mask the metal part. The SR Ivocron Opaquer is self-curing. Polymerization is completed after approx. 15 minutes.

Tip

Before proceeding, use a suitable instrument to check whether the material has cured completely.

Dentin and Incisal

Mix the desired SR Ivocron Dentin or Incisal powder in a rubber cup, cover it, and let it sit for approx. 3–4 minutes.

Tip

Mixing Dentin and Incisal materials in a 1:1 ratio results in improved shade adaption in interdental areas.

Pour the Ivocron mixture in the silicone key, while it is still flowable. Allow to stand briefly, and immediately place it in the pressure apparatus.

Polymerization

Polymerization is carried out in the pressure apparatus (e.g. Ivomat) for 15 minutes at 40–50 °C (104–122 °F) and 2-6 bar pressure (29–85 psi).

Finishing, polishing

Finish the restoration according to the standard procedure for PMMA materials, i.e. finishing with burs and rubber polishers, pre-polishing with goat's hair brushes and pumice, and polishing to a high gloss with cotton wool buffers and polishing medium (e.g. Ivoclar Vivadent Universal Polishing Paste).

Hot technique (layering technique)

When used together with Hot Liquid, SR Ivocron can be mixed to a smooth consistency. In this way, the veneer can be easily layered. The excellent material quality is achieved by polymerization in the pressure apparatus at a minimum temperature of 100 °C (212 °F). This technique is preferable for long-term temporaries and inexpensive crowns and bridges.

Fabricating the framework

Use metal to fabricate the palatal and incisal frameworks of anterior bridges or the occlusal area of posterior bridges. Apply mechanical retentions (e.g. SR Micro and SR Macro Retention Beads from Ivoclar Vivadent) on the veneering surface.

Isolation

Block out undercuts with pink wax. Then, immerse the model in water for 5 minutes. Apply two layers of Separating Fluid. Wait until the surface is no longer shiny before applying the second layer.

Bond

Use aluminium oxide to blast the metal parts of the bridge, which will be veneered. Since the opaquer does not chemically bond to metal, mechanical retentions have to be applied on the metal.

Tip

We recommend using SR Link for the chemical bond to metal. Blast framework with aluminium oxide and remove grinding dust. Wet it with SR Link and allow it to dry for 3 minutes.

Opaquer

Mix the desired opaquer with the SR Ivocron Opaquer Liquid, cover it, and let it sit for approx. 2–3 minutes. After that, mask the area to be veneered. The SR Ivocron Opaquer is self-curing. Polymerization is completed after approx. 15 minutes.

Tip

Before proceeding, use a suitable instrument to check whether the material has cured completely.

Dentin and Incisal

Mix the desired SR Ivocron Dentin or Incisal powder in a rubber cup, cover it, and let it sit for approx. 2–3 minutes. Replace the framework on the isolated model. Then, layer the veneer using a brush or other suitable instrument.

Polymerization

Polymerization is carried out in the pressure apparatus (e.g. Ivomat) for 30 minutes at 100 °C (212 °F) and 2–6 bar pressure (29–85 psi).

Tip

Off the model, polymerization can be carried out in the pressure apparatus (e.g. Ivomat) for 10 minutes at 120 °C (248 °F) and 2–6 bar pressure (29–85 psi). The stone would be destroyed at these temperatures.

Finishing, polishing

Finish the restoration according to the standard procedure for PMMA materials, i.e. finishing with burs and rubber polishers, pre-polishing with goat's hair brushes and pumice, and polishing to a high gloss with cotton wool buffers and polishing medium (e.g. Ivoclar Vivadent Universal Polishing Paste).

Pressing technique (flasking technique)

For the flasking technique, SR Ivocron is mixed with the Press Liquid. The material then demonstrates a smooth consistency and can be easily pressed.

Fabricating the framework

Use metal to fabricate the palatal and incisal frameworks of anterior bridges or the occlusal area of posterior bridges. Apply mechanical retentions (e.g. SR Micro and SR Macro Retention Beads from Ivoclar Vivadent) on the veneering surface.

Bond

Use aluminium oxide to blast the metal parts of the bridge, which will be veneered. Since the opaquer does not chemically bond to metal, mechanical retentions have to be applied on the metal.

Tip

We recommend using SR Link for the chemical bond to metal. Blast framework with aluminium oxide and remove grinding dust. Wet it with SR Link and allow it to dry for 3 minutes.

Opaque


Mix the desired opaquer with the SR Ivocron Opaquer Liquid, cover it, and let it sit for approx. 2–3 minutes. After that, mask the part to be veneered. The SR Ivocron Opaquer is self-curing. Polymerization is completed after approx. 15 minutes.

Tip

Before proceeding, use a suitable instrument to check whether the material has cured completely.

Modelling the crown shape/investing

Make an accurate model of the veneer with tooth-coloured wax. Invest the waxed-up model in the flask using dental stone.

Important

Observe an angle of 25–30° and make sure that the object offers proximal access.

Boiling out and isolation

Open the flask and boil out. Isolate the warm plaster halves with two layers of Separating Fluid. Wait until the surface is no longer shiny before applying the second layer.

Preheating the core

Preheat the core in boiling water to prepare it for the pressing procedure. In this process, the water must not touch the isolated surface.

Wetting the opaquer

Wet the exposed opaquer with monomer.

Mixing the Dentin material

Mix SR Ivocron Dentin material to a viscous consistency in a rubber cup, cover it, and let it sit for approx. 2–3 minutes. The material is ready for use as soon as it no longer sticks to the spatula.

Pressing the Dentin material

Apply Dentin material on the framework, slightly adapt it with a suitable instrument, and cover it with a plastic foil. Place the pre-heated second flask half on the first half. Apply 2 tons of pressure to the flask and allow to cool under pressure.

Reducing the incisal area

Open the flask, remove foil, and pre-heat the upper flask half in boiling water again. For the incisal area, bevel the dentin material with a sharp instrument.

Mixing the Incisal material

Mix SR Ivocron Incisal material in a rubber cup (the mixture should be less viscous than the dentin mixture) and let it sit for 3–4 minutes. Then, apply the material with a suitable instrument and spread it.

Pressing the Incisal material

Cover the material with a plastic foil. Place the preheated second flask half on the first half and close flask. Press immediately. Pressure is approx. 1.5 tons, i.e. somewhat lower than that needed for the dentin. Allow to cool under pressure.

Layering of the cervical area

Open the flask, remove foil, and preheat the upper flask half in water again. For the cervical area, bevel the dentin material with a sharp instrument.

Mixing the Cervical material

Mix SR Ivocron Cervical material in a rubber cup, cover it and let it sit for 2–3 minutes. The material is ready to use as soon as it no longer sticks to the spatula.

Applying the Cervical material, pressing

Apply the mixed Cervical material on the exposed cervical area and cover it with a plastic foil. Place the preheated second flask half on the first half, close the flask, and press immediately. Pressure is approx. 1 ton, i.e., somewhat lower than that needed for the incisal area. Allow to cool under pressure.

Polymerization

Secure the flask in the clamping frame.

– Single crowns

Place in cold water, heat, and boil for 30 minutes. Allow to cool slowly.

– Bridges

Preheat in water measuring 70° C (158 °F) for 30 minutes. Then, boil for 30 minutes, and allow to cool slowly.

Finishing and polishing

Finish the restoration according to the standard procedure for PMMA materials, i.e. finishing with burs and rubber polishers, prepolishing with goat's hair brushes and pumice, and polishing to a high gloss with cotton wool buffers and polishing medium (e.g. Ivoclar Vivadent Universal Polishing Paste).

Characterization

SR Ivocron Intensive powders are reactive polymers. They are used for strong characterizations or for shade modifications of Dentin, Cervical, or Incisal materials.

Depending on the type of the restorations, SR Ivocron Intensive materials are mixed with Cold, Hot, or Press Liquid. Make sure that the same liquid is used as for the other materials. Small quantities of dry material may also be applied and then wetted with monomer.

Summary of the most important data

Cold Technique

- *Mixing ratio by volume*
1 part polymer : 1 part monomer
- *Mixing ratio in g*
1 g polymer : 0.83 g monomer
- *Dough time*
3–4 minutes
- *Working time at 23 °C (73 °F)*
Approx. 8 minutes
- *Polymerization*
In the pressure apparatus at 2–6 bar pressure and 40–50 °C (104–122 °F) for 15 minutes

Hot Technique

- *Mixing ratio by volume*
1–1.5 parts polymer : 1.5 parts monomer
- *Mixing ratio in g*
1–1.5 g polymer : 1.25 g monomer
- *Dough time*
2–3 minutes
- *Working time at 23 °C (73 °F)*
Depending on the mixing ratio approx. 8–25 minutes
- *Polymerization with model*
In the pressure apparatus at 2–6 bar pressure and 100 °C (212 °F) for 25–30 minutes
- *Polymerization without model*
In the pressure apparatus at 2–6 bar pressure and 120 °C (248 °F) for 10 minutes

Pressing Technique

Dentin and Cervical materials

- *Mixing ratio by volume*
2.5 parts polymer : 1 part monomer
- *Mixing ratio in g*
2.5 g polymer : 0.83 g monomer
- *Dough time*
2–3 minutes
- *Working time at 23 °C (73 °F)*
Approx. 8–10 minutes

Incisal material

- *Mixing ratio by volume*
2 parts polymer : 1 part monomer
- *Mixing ratio in g*
2 g polymer : 0.83 g monomer
- *Dough time Incisal*
3–4 minutes
- *Working time at 23 °C (73 °F)*
Approx. 8–10 minutes

Polymerization of a single crown

In the flask in the clamping frame, starting in cold water until the water boils for 30 minutes and cooling slowly.

Polymerization of a bridge

In the flask in the clamping frame. Preheating at 70 °C (158 °F) for 30 minutes. Then boiling for 30 minutes and cooling slowly.

Storage

- Store out of the reach of children
- Observe storage information and date of expiration on the secondary packaging
- Do not use the products after the indicated date of expiration

Shades

With SR Ivocron, all 20 Chromascope shades can be reproduced. We offer 20 cervical materials that can be used in conjunction with the SR Vivodent PE shade guide. Moreover, 10 Intensive polymers are available for shade modifications.

Intensive

- | | |
|-------------------|----------------|
| 1 – clear | 6 – dark brown |
| 2 – white | 7 – pink |
| 4 – yellow-orange | 9 – bluish |
| 5 – light-brown | |

Date information prepared: 01/2008

Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.